


# Comfort Cell

## Technology

### What is "Comfort Cell" Technology?

Comfort Cell Technology is a block of softer density foam inserted into the seat under the base of the spine, surrounded by firmer density foam. This is the result of "in depth research" with Occupational Therapists and considerable development to find the perfect combination of foam density. The Comfort Cell cushions the base of the spine while the surrounding firmer density foam holds the correct sitting position preventing the pelvis from slipping forward into a "bad posture" position. The result is that a "good posture" sitting position is maintained effortlessly all day reducing back aches, headaches, shoulder and neck fatigue and promoting healthy well being and performance. The "Waterfall" seat front completes the therapeutic effect by facilitating blood circulation behind the thighs.



- 1 'Comfort Cell' foam supports the pelvis yet allows it to sink in. This foam also absorbs heat generated by sitting for long periods of time, this leads to less body perspiration.
- 2 Cold cured and injection moulded foam with lifetime elastic memory.
- 3 Water fall front facilitates blood circulation behind the thighs.
- 4 Firmer foam in front section helps prevent the pelvis from moving forward in the uncomfortable 'C' curve of poor posture.
- 5 Raised sides of seat allows the person seated to always be central on the chair and hence improves healthy well being and performance.
- 6 The seat pan is available in 4 separate sizes which enables us to meet everyone's needs. Extra large, Large, Medium, Small.